

# RADIOLOGY DEPARTMENT EQUIPMENT SPECIFIC TRAINING

(Digital Mammography Unit)



Hospital	
Equipment Name	
Operator	
Assessor	

KNOWLEDGE	Tick
Demonstrate how to switch the unit on and perform the warm-up procedure.	
Understand the Quality Assurance procedure and be able to act on the results.	
Be able to move the gantry in all directions.	
Configure the system for the required view and procedure.	
Understand the control console and all its functions	
Enter patient details / select patient from the work list.	
Understand all automatic modes and manual mode functions.	
Choose and select an exposure mode and target selection as appropriate .	
Understand how to use the compression of the unit, and the force applied.	
Have knowledge of the emergency compression release.	
Have knowledge of how to send images to PACS	
Find the Protocols, Local Rules and Exposure Charts.	
If a failed message comes up when exposing, describe your next action.	
Have knowledge of how to remove the room from action.	
Demonstrates an understanding of infection control.	
Find the emergency stop button & how to reset.	

I certify that the above person has been trained and is competent to use the equipment without supervision. I further certify that this list is not exhaustive and does not constitute comprehensive training.

Signed \_\_\_\_\_ Assessor \_\_\_\_\_ Date \_\_\_\_\_

I have received adequate training and I feel that I am competent to use this equipment without supervision.

Signed \_\_\_\_\_ Operator \_\_\_\_\_ Date \_\_\_\_\_